

# Pharmacy Accessory Label Printer (PALP)


A bar code and labeling solution for your ancillary packaging needs

- Quickly and efficiently prints accurate bar coded labels on ampules, vials, and syringes.
- Utilizes MPI's WinPakUD™ pharmacy software
- Many different bar code formats available
- User maintained product master file allows for quick recall of standard labels
- Security options allow control of data files and labeling options
- Small footprint
- Low maintenance
- Four (4) different label sizes available.




Pharmacy Accessory Label Printer


## Pharmacy Accessory Label Printer Sample Thermal Labels



Item #DTLP-L1 Direct Thermal Labels (1380 Labels/Roll)  
1"x3" Size (3 Rolls per Case)



Item #DTLP-L2 Direct Thermal Labels (2400 Labels/Roll)  
1/2" x 3" Size (3 Rolls per Case)



Item #DTLP-L5 Direct Thermal Labels (2500 Labels/Roll)  
1/4" x 1-1/4" Size (4 Rolls per Case)



Item #DTLP-L4 Direct Thermal Labels (2750 Labels/Roll)  
1/2" x 1-1/2" Size (3 Rolls per Case)

### Machine Specifications

Width: 7.6"/103mm
Depth: 8.5"/216mm
Height: 6.6"/168mm
Weight: 3.2 lbs/1.5kg (weight of printer only; no media or options)