



## PHARMACY ACCESSORY LABEL PRINTER (PALP) – FACT SHEET

The **Pharmacy Accessory Label Printer** quickly and efficiently prints accurate bar coded labels on ampoules, vials and syringes. The PALP is comprised of a Zebra clamshell printer coupled with MPI's WinPak™ Barcode Labeling Software and is currently offered in two (2) different models to best meet the specific needs of each facility:

- **Direct Thermal Label Printer** – Provides the capability to print six (6) different size white labels, including white flag labels.
- **Thermal Transfer Label Printer** – Provides the capability to print clear stem flag labels only.

Both printers provide the following capabilities:

- Multiple bar code formats
- User-maintained product master file for quick recall of standard labels
- Security options to all control of data files and labeling options.
- Small footprint
- Includes MPI's exclusive WinPak™ Barcode Labeling Windows™ based software package.

### SPECIFICATIONS

<b>Width:</b>	7.6"/103mm
<b>Depth:</b>	8.5"/216mm
<b>Height:</b>	6.6"/168mm
<b>Weight:</b>	3.2 lbs./1.5kg (weight of printer only; no media or options)
<b>Hardware Reqs.:</b>	Computer with Windows® 98 or higher, 30 MB hard drive space, 9-pin serial port, desktop or network printer for reports, and optional Internet access for future software updates.
<b>Software Reqs.:</b>	WinPak™ Barcode Labeling Software (provided with packaging system)
<b>Pre-Installation:</b>	The PALP requires one (1) 120V outlet. Allow additional room and electrical outlets for your computer.

### PRICING

Item #	Description
DTLP-P1	<b>Direct Thermal Label Printer – (WHITE LABELS)</b>
TTLP-P1	<b>Thermal Transfer Label Printer – (CLEAR STEM FLAG LABELS ONLY)</b>

### DISPOSABLE MATERIALS

The **Direct Thermal Label Printer** requires one type of labels (see selection below) for standard operation. The minimal amount of materials required to operate the packaging system must be ordered at the time the equipment is purchased to ensure that proper testing can be performed and full functionality achieved at the time of installation and training.

Item #	Description	Labels/Roll	Qty/Case
DTLP-L1	1" x 3" Direct Thermal Labels – WHITE	1380	3 rolls
DTLP-L2	½" x 3" Direct Thermal Labels – WHITE	2400	3 rolls
DTLP-L3	¼" x 1" Direct Thermal Labels – WHITE	2500	4 rolls
DTLP-L4	½" x 1-½" Direct Thermal Labels – WHITE	2750	3 rolls
DTLP-L5	¼" x 1-¼" Direct Thermal Labels – WHITE	2500	4 rolls
DTFLAG-1	1-½" x 1 Direct Thermal Flag Labels – WHITE	1000	3 rolls

The **Thermal Transfer Label Printer** requires labels and thermal transfer ribbon for standard operation. Materials are currently being offered in a standard combo pack including both the labels and ribbon. Additional thermal transfer ribbons can be ordered separately. The minimal amount of materials required to operate the packaging system must be ordered at the time the equipment is purchased to ensure that proper testing can be performed and full functionality achieved at the time of installation and training.

Item #	Description	Labels/Roll	Qty/Case
TTFLAG1	Clear Stem Flag Label Combo Pack includes: <ul style="list-style-type: none"> <li>• Clear Stem Flag Labels</li> <li>• Thermal Transfer Ribbon</li> </ul>	1000 4.3" x 2904"	3 rolls 1 roll
TTLP-TR104	Thermal Transfer Ribbon	4.3" x 2904"	1 roll