



# MPI Auto-Wrap® Syringe Labeling System – Disposable Materials Required for Labeling System Operation

*MPI's Auto-Wrap® Syringe Labeling System requires any style Auto-Wrap® Label (shown below) and transfer ribbon to print, wrap and affix labels to a syringe.*

*Three (3) different styles of labels and four (4) different colored labels are available to meet the specific syringe labeling needs of your facility. Label print area ranges from 1.0" to 1.5".*




## FLAGGED LABELS



The top clear portion of the label is wrapped around the syringe with semi-permanent adhesive. The white/printed portion extends and hangs from the syringe. There is no adhesive behind the white area.

## TACKED LABEL



The top clear portion of the label is wrapped around the syringe with semi-permanent adhesive.

The white/printed portion has a clear adhesive area which will keep the label tacked to the syringe. This will allow for better storage and for the label to be untacked from the syringe for viewing.

**Please see reverse side for a list of the MPI disposable materials available for use with this labeling system.**

*Approximate label sizes shown. Materials are not shown at actual size.*

020117



**MEDICAL PACKAGING INC.**

8 Kings Court | Flemington | New Jersey | 08822 | USA  
Toll-Free.800-257-5282 | Phone.609.466.8991 | Fax.609.466.3775 | [www.medpak.com](http://www.medpak.com)

Auto-Wrap® is a trademark of Medical Packaging Inc.

**MPI MATERIALS AVAILABLE FOR USE WITH  
AUTO-WRAP® SYRINGE LABELING SYSTEM**

| <b>GTIN Number</b> | <b>Item #</b> | <b>Item Description</b>  |
|--------------------|---------------|--|
| 00818259011056     | SL-11510      | <b>Auto-Wrap®</b> Syringe Label ( <b>Flag</b> ), 3" x 1.5" ( <b>1.0"</b> print area) - 6,000 (6x1000)    |
| 00818259011063     | SL-11515      | <b>Auto-Wrap®</b> Syringe Label ( <b>Flag</b> ), 3" x 1.5" ( <b>1.5"</b> print area) - 6,000 (6x1000)    |
| 00818259011070     | SL-21512      | <b>Auto-Wrap®</b> Syringe Label ( <b>Tacked</b> ), 3" x 1.5" ( <b>1.25"</b> print area) - 6,000 (6x1000) |
| 00818259011155     | SL-31515G     | <b>Auto-Wrap®</b> Syringe Label (Flag), Green, 3" x 1.5" (1.5" print area) - 6,000 (6x1000)              |
| 00818259011186     | SL-31515O     | <b>Auto-Wrap®</b> Syringe Label (Flag), Orange, 3" x 1.5" (1.5" print area) - 6,000 (6x1000)             |
| 00818259011179     | SL-31515R     | <b>Auto-Wrap®</b> Syringe Label (Flag), Red, 3" x 1.5" (1.5" print area) - 6,000 (6x1000)                |
| 00818259011162     | SL-31515Y     | <b>Auto-Wrap®</b> Syringe Label (Flag), Yellow, 3" x 1.5" (1.5" print area) - 6,000 (6x1000)             |
| 00818259010813     | T002P         | Ribbon, Trnfr SL 2.09" - 7,080LFT (6x1181')  |