

# PHARMACY ACCESSORY LABEL PRINTER (PALP)

Manufactured by Zebra Technologies



DECREASE PACKAGING COSTS

INCREASE OPERATIONAL EFFICIENCY

REDUCE MEDICATION ERRORS

## Pharmacy Accessory Label Printer (PALP)

is an easy-to-use labeling system for printing labels for parenteral medications with MPI's exclusive Pak-EDGE™ UD Barcode Labeling Software.



Zebra GK420 Series

MPI's **Pharmacy Accessory Label Printer (PALP)** provide a low cost solution to quickly and efficiently print accurate bar-coded labels for placement on ampoules, vials and syringes. MPI's PALP solutions, comprised of a Zebra GK 420d printer and MPI's exclusive **Pak-EDGE™ UD Barcode Labeling Software**, are offered with either direct thermal or thermal transfer technology to pharmacy:

- **Direct Thermal Label Printer** – Provides the capability to print standard, circular, and white flag labels. Yellow, red and orange flag labels are also available.
- **Thermal Transfer Label Printer** – Provides the capability to print clear stem flag labels only; two flag label options available.

MPI's **Pharmacy Accessory Label Printer (PALP)** is powered by MPI's exclusive **Pak-EDGE™ UD Barcode Labeling Software**, which is capable of generating and printing all linear, 2D and GS1 barcodes and offers advanced barcoding capabilities and configurations, which typically includes product NDC code, local expiration date and date packaged. In the event of a product recall, Pak-EDGE™ can easily trace the packaged product so the pharmacist can take immediate corrective action to remove the product from inventory.

Labeling parenteral medications with MPI's **Pharmacy Accessory Label Printer (PALP)** produces a legible, ready to scan label containing detailed medication and barcode information assuring that the correctly labeled medication is delivered to the bedside each and every time, ensuring patient safety.

### KEY FEATURES

- Multiple bar code formats including all 2D and linear formats and Advanced Barcode Builder capability.
- User-maintained product master file allows for quick recall of standard labels.
- Enhanced security options allow control of data files and labeling options.
- 203 dpi/8 dots per mm resolution.
- Minimal Cost Per Dose (CPD)
- Compact footprint

### TECHNICAL SPECIFICATIONS

Dimensions:	8.25"L x 6.75"W x 6.0"H
Weight:	3.0 lbs.
Safety:	CE, UL, C
Power:	100-240VAC, 50/60Hz

All MPI packaging systems are classified as US FDA 510(k) exempt and require the use of MPI authorized disposable materials.



To request information, contact:  
[sales@medpak.com](mailto:sales@medpak.com)

 **Medical  
Packaging Inc.**

8 Kings Court | Flemington | New Jersey | 08822 | USA  
Toll-Free.800-257-5282 | Phone.609.466.8991 | Fax.609.466.3775 | [www.medpak.com](http://www.medpak.com)



02/01/17

Pak-EDGE is a trademark of Medical Packaging Inc.