



MPI Auto-Wrap® Syringe Labeling System – Disposable Materials Required for Labeling System Operation

Medical Packaging Inc., LLC (MPI)'s *Auto-Wrap® Syringe Labeling System* requires any style *Auto-Wrap® Label* (shown below) and transfer ribbon to print, wrap and affix labels to a syringe.

Three (3) different styles of labels and five (5) different colored labels are available to meet the specific syringe labeling needs of your facility. Label print area ranges from 1.0" to 1.5".




FLAGGED LABELS



The top clear portion of the label is wrapped around the syringe with semi-permanent adhesive. The white/printed portion extends and hangs from the syringe and the area is deadened to reduce tackiness.

TACKED LABEL



The top clear portion of the label is wrapped around the syringe with semi-permanent adhesive.

The white/printed portion has a clear adhesive area which will keep the label tacked to the syringe. This will allow for better storage and for the label to be untacked from the syringe for viewing.

Please see reverse side for a list of the MPI disposable materials available for use with this labeling system.

Approximate label sizes shown. Materials are not shown at actual size.



8 Kings Court | Flemington | New Jersey | 08822 | USA
Toll-Free.800.257.5282 | Phone.609.466.8991 | Fax.609.466.3775 | www.medpak.com

Auto-Wrap® is a trademark of Medical Packaging Inc., LLC.

MEDICAL PACKAGING INC., LLC (MPI)
DISPOSABLE MATERIAL AVAILABLE FOR USE WITH
MPI's AUTO-WRAP® SYRINGE LABELING SYSTEM

GTIN Number	Item Number	Item Description
00818259011056	SL-11510	Auto-Wrap® Syringe Label (Flag), 3"x 1.5" (1.0" print area) - 6,000 (6x1000)
00818259011063	SL-11515	Auto-Wrap® Syringe Label (Flag), 3"x 1.5" (1.5" print area) - 6,000 (6x1000)
00818259011070	SL-21512	Auto-Wrap® Syringe Label (Tacked), 3"x 1.5" (1.25" print area) - 6,000 (6x1000)
00818259011155	SL-31515G	Auto-Wrap® Syringe Label (Flag), Green, 3"x 1.5" (1.5" print area) - 6,000 (6x1000)
00818259011186	SL-31515O	Auto-Wrap® Syringe Label (Flag), Orange, 3"x 1.5" (1.5" print area) - 6,000 (6x1000)
00818259011179	SL-31515R	Auto-Wrap® Syringe Label (Flag), Red, 3"x 1.5" (1.5" print area) - 6,000 (6x1000)
00818259011216	SL-31515V	Auto-Wrap® Syringe Label (Flag), Violet, 3"x 1.5" (1.5" print area) - 6,000 (6x1000)
00818259011162	SL-31515Y	Auto-Wrap® Syringe Label (Flag), Yellow, 3"x 1.5" (1.5" print area) - 6,000 (6x1000)
00818259010813	T002P	Transfer Ribbon, 2.09" - 7,086LFT (6x1181')